




Status of MXCuBE Beamline Control at BESSY II

Michael Hellmig,
on behalf of the HZB-MX group

MXCuBE/ISPyB Joint Meeting, 29.10.-31.10.2019,
BESSY II @HZB, Berlin

Synchrotron sources at BESSY II




7 Tesla wavelength shifter and MX Beamlines

BESSY II ring parameters:


Electron Energy:	1.7 GeV
Electron Current:	280 mA
Circumference:	240 m
Straight sections:	16
Beamlines:	~50

main building „Helmholtz-Zentrum Berlin“



HZB-MX BEAMLINES OVERVIEW


U. Mueller *et al.* (2012).
J. Synchrotron Rad. **19**, 442-449.
U. Mueller *et al.* (2016).
Eur. Phys. J. Plus **130**, 141-150.



optics hutch


experimental hutches

user control room



More than 80 groups
from 15 countries

source: biosync.sbkb.org



*) 2018-depositions will be complete 12/2019

MX experimental floor at BESSY II

BL 14.1 MAD

- MD2 with MK3
- Pilatus2 6M 12 Hz
- CATS: 90 SPINE samples
- MXCuBE 2.2 Qt4




- standard user operation schedule:
24/5 (Tuesday to Saturday)



BL 14.3 13.8 keV

- MD2S with MK3
- Rayonix MX225
- HClab & REX nozzle changer
- MXCuBE 2.2 Qt4



user operation

BL 14.2 MAD

- Nanodiff goniometer
- Pilatus3 2M
- GROB: 294 SPINE & Unipuck samples
- MXCuBE 2.2 Qt4




back in operation




BL 14.1 MAD

- detector upgrade planned for 2020
 - existing detector installed 01/2013
 - manufacturer support ended
 - no issues with detector head
 - detector control computer: Ubuntu 10.04 LTS
 - replacement of Pilatus2 6M with current model Pilatus3 6M (S or X)
 - PPU for fast online data analysis
- re-use existing detector at BL14.3
 - replace Rayonix CCD detector
 - enable shutterless data acquisition
 - ~3 times faster data collection
- CATS Uni-puck upgrade under evaluation
 - SPINE and Uni-puck support
 - shorter exchange time with Uni-puck double gripper



- upgrade of beam-shaping devices in progress
 - unreliable operation
 - slow
 - on-call service intervention



short travel range




insufficient isolation against LN2

BL 14.2 MAD

- new setup by Smaract
 - larger travel range
 - transparent control-system integration with Tango DS
 - to be installed 12/2019




MXCUBE-HCLab integration



BL 14.3 13.8 keV

- integration in progress
 - HCLabGroupQueueEntry
 - HCLabQueueEntry
 - HCLabDataCollectionQueueEntry



Oleg Kornelsen
Bachelor student

- TO-DO:
 - HCLab status widget
 - online analysis: DISTL.Spotfinder/Dozor

MXCuBE 2.2 Qt4 experiment-control software at all three HZB-MX beamlines

Diffraction meters:

Arinax MD2
Arinax MD2S
DESY Nanodiff

Sample-transfer robots:

Irelec CATS
NatX-ray GROB

Detectors:

Dectris Pilatus2 6M
Dectris Pilatus3 2M
Rayonix MX225

Auxiliary devices:


Wago I/O controller
Amptek X-123SDD

Control systems:

EPICS, Tango,
Exporter, SPEC

MXCuBE: status and future plans

- current status:
 - MXCuBE 2.2 Qt4 running on all three HZB-MX beamlines
- current issues:
 - resources mostly spend on testing, maintenance, trouble shooting
 - commissioning/testing time at the beamline strongly limited
 - outdated MXCuBE control computer at the beamline
 - Debian 7, **32-bit**
- short- and mid-term plans:
 - modernization of MXCuBE control computer
 - Debian 9/10, 64-bit
 - Anaconda/Miniconda development environment
 - offline MXCuBE update:
 - unified HZB-specific HardwareObjects „branch“
 - HardwareRepository → branch 2.3.0 → master
 - MXCuBE2 → GitHub master branch
 - complete integration of HCLab into MXCuBE software setup
 - HCLab brick
 - analysis



FragMAX webapp: collaboration MAX IV and HZB:

**identification
of hit fragments**

**auto-refinement
Pipeline
and PanDDA**

The diagram shows a 3x3 grid of panels (A-I) illustrating different refinement techniques: A: MR (3.0 or 2.5 Å), B: Simulated annealing, C: Standard refinement, G: Anisotropic B factors (water), H: Hydrogen addition, I: Additional water placement.

XDSAPP


The image shows the XDSAPP software interface with various plots and data processing options.

**automated data
processing**

Gustavo Lima (MAX IV), Jan Wollenhaupt (HZB)



XtalTool:

- All-in-One sample holder
- *in situ* method
- Crystallize on sample holder
- Use as a lid
- Data collection at beamline with magnetic base
- Developed by [C. Feiler](#), patent
- Distributed by Jena Bioscience



X-ray transparent Kapton
film with 5 μm pores

Transparent
removable COC film



BESSY-MX team

Manfred Weiss

Christian Feiler

Ronald Förster

Martin Gerlach

Christine Gless

Thomas Hauß

Huiling He

Michael Hellmig

Alexandra Kastner

Frank Lennartz

Michael Steffien

Helena Taberman

Jan Wollenhaupt



The MXCuBE collaboration



Industrial partners:



Thank you for your attention.

Questions?